

Make a step change

If a Stannah stairlift could make life easier, why wait?

Life is better when you have the energy to do the things you want to do, like enjoying a visit from the grandchildren or entertaining your family and friends.

Our experience has taught us that the sooner a stairlift is installed, the sooner you will be able to start enjoying the things you love doing.

Just think how a Stannah stairlift could make your day-to-day life easier.

A stairlift that's right for you and your home.

We know that not all homes are the same. If you think your stairs are unsuitable for a stairlift, we think you might be pleasantly surprised by our slimline and stylish designs.

- Our slimline rails and fold-away chairs make the best use of available space and leave enough room for others to use the staircase.
- All our rails are tailor-made just for you, so they fit your staircase perfectly.
- Contrary to popular belief our rails are fitted to the staircase and not the wall ensuring minimal mess and bother.
- Installation is quick and easy.



A gentle touch of the controls is all it takes to climb your stairs

A company you can trust

We are a British, family-run and owned company and have been for five generations.

The current managing directors are the great, great grandchildren of Joseph Stannah, who started the company in the 1860s. His engineering principles and strong family values run as true today as they always have.


Now, over 150 years later, we have manufactured more than 700,000 stairlifts and installed products in over 40 countries.

Starla: Stylish, sleek and very comfortable

The elegant Starla is crafted for choice and ease-of-use in every way. The chair folds away for a compact footprint an automatic footrest fold button so you don't need to bend down.

The Starla comes with a choice of chair upholsteries and combinations so you can make it your own.

Tailor it to your tastes and home by choosing a fabric, with or without a light or dark wood finish or select the full vinyl upholstery in a range of sophisticated colours.



The Starla's slimline design allows others to still use the stairs

Starla key features

1. Seat height

The seat height can be set to the most comfortable setting for you when the stairlift is installed, allowing you to get on and off the chair easily.

2. Wall controls

Ideal for when you need to summon your stairlift from another floor, or perhaps send it up or down to someone else who wants to use it.

3. Footrest

For added peace of mind, the footrest comes with an anti-slip material and a soft fold mechanism.

4. Arm isolation

For added safety, the stairlift won't be able to move until the armrest closest to the ground floor is fully lowered.



5. Easy-to-use seatbelt

To offer you even more security, we have specifically designed a seatbelt that can be easily fastened with one hand.

6. Safety edges

In the event that something obstructs the stairlift, sensors on the edge of the carriage and footrest will automatically bring the stairlift to a safe stop.

7. Easy-to-use controls

Effortlessly ride the Starla using our simple, intuitive controls. Choose between a paddle or rudder style which work with the lightest of touches, and can be operated using your hand, arm or wrist.

Starla key benefits


Effortless and simple. This sums up beautifully just how much thought and care we have designed into each and every function of our Starla stairlift. From the safe and easy-to-use seatbelts, to stairlifts that take up minimal space on your staircase, we have made your Starla stairlift simple, safe and reliable.



Seatbelt can be fastened with one hand



Use the wall control to call your Starla



Sensors detect any obstructions on your stairs

"We are delighted with the stairlift. My wife enjoys the ride. It gives her independence. She finds it very easy to use"

Mr A, Cirencester



The easy way to climb your stairs

Everything about a Stannah stairlift is made to make managing your stairs as easy, comfortable and safe as it can be, so your life changes for the better the moment you begin to use it.

From the simple controls at your fingertips to reliable operation and a flexible and compact design, we have done all the hard work so you don't have to. Just settle into the comfortable seat and sit back, as you travel up or down your stairs with ease.

We really have taken care of every detail to ensure your Starla stairlift is as effortless to use as it is easy to operate.

Step 1: Simply pull down the arms and then the seat base and the footrest will safely unfold into position.

Step 2: When the seat base and footrest are down, sit on your stairlift and make yourself comfortable in the

seat. Then pull across and fasten the easy reach seatbelt and you're now ready to go.

Step 3: A gentle movement of the directional control is all it takes for you to travel effortlessly and safely up or down your stairs.


Step 4: When you have arrived at the top of your stairs, you can then use your

swivel lever or optional powered swivel to rotate you onto your landing ready to get off your Starla in a safe position.

Step 5: When you are not using your stairlift, it can remain in position at the top of your stairs to act as a safety barrier, or, it can be folded away to allow space for other stair users.

"Now we have our Stannah stairlift, we no longer have to face an unwanted move into single storey accommodation"

Mr & Mrs H, Worcestershire











Tailored to you

The Starla is our most customisable stairlift. In addition to the upholstery options, you can also choose from a range of other features that will ensure your Starla is just right for you.


Seatbelt choices







Diagonal seatbelt



Five point harness

Control choices and control position



Light-touch control

Left hand Right hand



Rudder style control

Folding footrest



Powered footrest button

Swivel seat options



Manual swivel seat



Powered swivel seat



Two-way powered swivel seat

Starla upholstery and rail colour choices

You can truly make your stairlift unique by selecting one of our many upholstery options. Along with our standard offering (shown below), you can choose from our wide range of collections, including the Agua Fusion Collection with its contemporary diamond pattern, our luxurious velvety fabric range and our stylish yet practical vinyl/fabric combinations.

Fully upholstered



Green Woven



Red Woven





Brown Vinvl



Stone Woven

Light Beige Vinyl

Light wood trim



Green Woven



Red Woven



Stone Woven










Dark wood trim













Silver



guidance only. Although we have tried to display the colours accurately, we cannot guarantee that the colours in the pictures accurately reflect the colour of the actual stairlifts. Your stairlift may be slightly different from these pictures.

The pictures within our brochures are for

Green Woven Red Woven

Stone Woven

Brown Vinyl

Light Beige Vinyl

A professional installation

When you have chosen the Starla that suits you best, a Stannah trained engineer will visit your home and professionally install your new stairlift in a matter of hours, leaving your home just as they found it.





Our company policy is one of continuing product improvement. We therefore reserve the right to alter the specifications and to change and improve design at any time. Full operating instructions are supplied for all products. Stannah is a registered trademark. The CE mark indicates that this product complies with the applicable European directives which relate to health, safety, environment and consumer protection. Design, manufacture and installation: Stannah Stairliffs Ltd has been accredited to ISO9001:2008 since 2009, this is for design, development, production and distribution of stair lifts. Alterations in models, standard and optional equipment as described in the text and illustrations, may occur. For precise information please contact Stannah.

© Stannah Stairliffs Ltd 2020. Not to be reproduced wholly or in part without the written permission of Stannah Stairliffs Ltd Issue 1 / January 2020

